




# Fast Photo-Curable Silicones for Additive Manufacturing and Organ-on-Chip Technologies

Julien Gautrot








Rubber in Engineering Group  
*4<sup>th</sup> December 2020*




# Engineering the Cell Microenvironment at Multiple Scales





Kong et al. *ACS Nano* 2018  
Kong et al. *Nano Letters* 2018  
Trappmann, Gautrot et al. *Nat. Mater.* 2012  
Costa et al. *Acta Biomater.* 2014  
Gautrot et al. *Nano Letters* 2014



100 nm




1 μm





Colak et al. *Biomacromolecules* 2018  
Di Cio et al. *Acta Biomater.* 2016  
Di Cio et al. *Acta Biomater.* 2017  
Tan et al. *Integ. Biol.* 2013  
Connelly et al. *Nat. Cell Biol.* 2010  
Gautrot et al. *Biomaterials* 2012

# Cells in Culture Do not Look Like Cells *in Vivo*

## 2D In Vitro Culture



## In Vivo





Niessen et al. *J. Cell Sci* (1996), *109*, 1695.



Herle et al. *Development* (1991), *112*, 193.

- Cells and tissues *in vivo* have reproducible shapes, size and geometries.
- In 2D cultures, cell shape is unconstrained.


## Skin - Epidermis



## Intestine








- Stem cells reside in well defined locations.
- The architecture of the tissue and that of the “niche” are important to ensure proper functioning of stem cells and tissue homeostasis.





# Recreating Higher Level Structure and Function




*Science* **328**, 1662 (2010)

- Recreates structure and captures biophysics.
- Probes observed nanotoxicological response to nanoparticles (immune response), in particular in biophysical context (mimicking breathing).





# Higher Degree of Structure and Function in Tissues






# An Engineering Approach






# Multi-Scale Engineering of Compartmentalised Microvascularised Tissues


ECM Biochemistry  
Physics/Mechanics  
In situ gellation




Biomimetic basement membrane






Compartmentalisation  
Microfluidics  
3D Printed




Bonding  
Long-term structure



Stabilised co-culture




Biomechanics  
Pneumatic





Perfusable vasculature



# Microfabrication and Prototyping of Photocurable Silicones








**Compartmentalisation  
Microfluidics  
3D Printed**



**Biomechanics  
Pneumatic**



# 3D Ink Design – Ultra-Fast Cure Chemistry




- Ultra-fast cure chemistry.
- Cures in air without any oil residue.
- No toxic catalyst.





# Graphene-Based Conductive Silicones

## Excellent Thermal Stability




- Thermal stability up to 400°C
- Curing of composites with 3 wt% graphene oxide under 1s (250 μm samples).
- Achieve high conductivities even at 0.5 wt% after Gox thermal conversion.

## Fast Gelation of Opaque Composites



## Flowable, Printable Conductive Silicone Composites



# 3D Printed Silicones – Flexible 3D Design for Organ-on-Chip Applications




- Flexibility of design and automation.
- True 3D microstructuring.
- Limited resolution ( $> 100$  microns).





# 3D Ink Design - Formulation

## Tuning Ink Mechanics




## Rapid Curing of Filled Formulations




- Design of rheological properties of uncured formulations and mechanical properties of cured resins independently.
- Retention of fast curing even with opaque samples.
- Thixotropic properties suitable for extrusion-based printing.

## Thixotropic Properties



# Chip Design and Printing





- More compact and symmetrical design to avoid uncontrolled nutrient/growth factor gradients.
- Introduce pneumatic chamber in upper compartment.




# Pneumatic Chamber for Biomechanical Actuation


- Simple pneumatic chamber connected to software controlled pump.
- *In silico* prediction of the deformability of the compartment.



## Predicting Deformation




## Cyclic Deformation





# Mechanically Integrated Biomimetic Hydrogels

ECM Biochemistry  
Physics/Mechanics  
In situ gellation






Bonding  
Long-term structure






# Mimicking the Cell Microenvironment




- Capture biochemical composition (cell adhesion, degradation, matrix deposition).
- Stiffness.
- Viscoelasticity.
- Porosity / morphology.

# Peptide-Based Cell Degradable Hydrogels



Y. You, K. Suzuki, J. Gautrot et al. *Biomaterials* 2020, 120356




- Fast curing in physiological conditions.
- Control of mechanical properties and degradability.







# Regulation of Cell Spreading and Secretory Phenotype

Day 14




Y. You, K. Suzuki, J. Gautrot et al. *Biomaterials* 2020, 120356

POx2-75






- Increased degradability correlates with cell spreading.
- Restriction of cell spreading correlates with increased growth factor secretion.
- Matrix engineering regulates pro-angiogenic phenotype.

# Hydrogel Bonding in Thiol-Ene Hybrids



- Thiol-ene based hydrogels and silicones are chemically compatible.
- Improved adhesion compared to Sylgard PDMS.





Stretching rate : 1mm/min




# Engineering of Stable Perfusable Microvasculatures in Microfluidic Chips



Stabilised co-culture




Perfusable vasculature



# MSC Co-culture Promotes the Formation of a Mature Microvasculature



F-Actin, CD31, DAPI




- Formation of perfusable microvascular network in chips.
- Interfacing with advanced multi-cellular in vitro models for safety/efficacy testing.

*Dibble, Luo, Jones*






# Pericyte Co-cultures for Stabilisation of Microvasculatures




Attwell et al. Journal of Cerebral Blood Flow & Metabolism (2016).

## Impact on Microvasculature Structure






# Long-Term Stability of Vascular Networks in Microfluidic Chips

HUVECs + Pericytes



Fibronectin  
CD31  
F-actin




Fibronectin

- Pericyte co-culture allows to stabilise microvascular networks for >3-4 weeks.
- Stable in multiple types of culture medium, even upon serum starvation.
- Compatible with implantation of more complex models.



## Conclusions

- Thiol-ene PDMS display excellent properties for 3D printing in ambient conditions.
- Controlled mechanics and interfacing with hydrogels.
- Enables the fabrication of 3D chips for biomechanical actuation and the embedding of large complex multi-cellular tissue models.





# Thank You



## Collaborators

John Connelly (QMUL)  
Fran Balkwill (QMUL)  
Wen Wang (QMUL)  
Adrian Biddle (QMUL)  
Madeline Lancaster (LMB, Cambridge)  
Pavel Novak (QMUL)  
Madeleine Ramstedt (Umea Universtiy)  
Armando del Rio (Imperial College)  
Ali Zarbakhsh (QMUL)  
Ken Suzuki (QMUL)  
Helena Azevedo (QMUL)

Duncan Sutherland (Aarhus)  
Jenny Malmstrom (Aarhus)  
Christoph Ballestrem (Manchester)  
Thomas Iskratsch (QMUL)  
Mariah Sundh (Aarhus)  
Fiona Watt (KCL)  
Wilhelm Huck (Nijmegen)  
Arnoud Sonnenberg (NKI)  
Thereza Soares (Recife)

Burcu Colak  
Stefania di Cio

Pei Tang

Dexu Kong

Khai Nguyen

Alex Moothedathu

Alice Dowden

Fengqin Qu

Mahentha Krishnamoorthy

Cardee Alexis

Lan Chang

Elijah Mojares

Danyang Li

William Megone

Matt Dibble

Piaopiao Luo

Ed Cozens

Yaqi You

Linke Wu

Christian Jones

Lihui Peng

Fauzia Quadir

Minerva Bosch

Alexandra Chrysanthou

Nam Nguyen

Chiara da San Martino




STINT

Stiftelsen för internationalisering av  
högre utbildning och forskning

The Swedish Foundation for International  
Cooperation in Research and Higher Education



National Centre  
for the Replacement  
Refinement & Reduction  
of Animals in Research



[biointerfaces.qmul.ac.uk](http://biointerfaces.qmul.ac.uk)

Positions available [j.gautrot@qmul.ac.uk](mailto:j.gautrot@qmul.ac.uk)